



# Instructions For Use

---













## **Tumour Protein p53 gene (*TP53*) Library Prep Kit for Illumina®**

NGS Library Prep by Reverse Complement PCR

---

*Version:* 2  
*REF:* IFU-TP53  
*Revision Date:* 2022-10-05

## Description of the symbols used on the labels

Symbol	Description
	Manufacturer
	Use-by date
	Lot number
	Reference number
	Upper limit of temperature for storage
	Temperature limit for storage
	Consult the instructions for use
	Contains sufficient for <n> tests
	CE-mark
	<i>In vitro</i> diagnostic medical device
	The box contains 2 index plates
	Matrix code containing the reference number, lot number and use-by date

## Product and Company Information

Product name: Tumour Protein p53 gene (*TP53*) Library Prep Kit for Illumina®  
REF: RC-TP53096-I  
Product use: For *In Vitro* Diagnostic use  
Company: NimaGen B.V.  
Hogelandseweg 88  
6545 AB Nijmegen  
The Netherlands  
Telephone: +31 (0)24 820 02 41  
Email: info@nimagen.com

## Product Use

The intended purpose of this assay is NGS library preparation for Illumina® sequencing of the human *TP53* gene coding exons, including a minimum of 20 bases upstream and downstream of each exon; to aid in diagnosis and or prognostication of cancer by detecting acquired mutations within the *TP53* gene; to aid in diagnosis of congenital predisposition to cancer by detection of germline mutations in the *TP53* gene (Li-Fraumeni syndrome). The lower detection limit of this assay is 5% variant allele frequency.

The assay is designed for use with DNA extracted from peripheral blood or tumour tissue. Specimens should have an absorbance ratio (260/280) of ~1.8, regardless of the extraction method used.

The assay provides reagents for Multiplex Amplicon-based NGS library preparation and is for professional use only.

The kit is based on the patented Reverse Complement PCR technology, providing a safe, robust and simple workflow, combining multiplex amplification in a single reaction, decreasing the risk of PCR contamination and sample swapping.

## Kit Content: RC-TP53096-I (-20 °C storage)

Description	Content	
RC-PCR Probe Panel A (REF: PM-TP53-A)	Tube 24 µL	●
RC-PCR Probe Panel B (REF: PM-TP53-B)	Tube 24 µL	●
2x Master Mix (Hot Start HiFi) (REF: MMHS096)	2 Tubes 1.15 mL	○
Probe Dilution Buffer (REF: RC-PDB)	Tube 500 µL	●

## To be ordered separately: IDX096-U01D-I or IDX096-U02D-I (max 25 °C storage)

Description	Content
IDX* Primer Plates, dehydrated, coloured, 96 Unique Dual, 10 bp Indexes for Illumina®	2x 96-well plate, with Alphanumeric Coding
Caps	2x 12 Domed 8-Cap Strips

*\*IDX plates are semi-skirted, ABI style PCR Plates, breakable per 8 tubes, containing dehydrated and stabilized Unique Dual Index primer pairs in each well and a red tracking dye, ready to use.*

The reagents do not contain any tissues or substances from human or animal origin.

## Needed, but not included

Description	Vendor
Adjustable Pipette Set (P2, P10, P20, P100, P200, P1000)	Multiple Vendors
TapeStation, Bioanalyzer System, or equivalent, incl. consumables, or optional: agarose gel system	Agilent® or other
Ethanol absolute, mol. biol. grade	Multiple Vendors
AMPureXP™ Bead Solution (or similar)	Beckman Coulter
General plasticware, DNase free (1.5 mL tubes, pipette tips w/filter)	Multiple Vendors
Ice or tabletop cooling block	Multiple Vendors
Mini Spinner for 1.5 mL tubes and 8-well PCR strips	Multiple Vendors
Magnetic Stand for 1.5 mL Eppendorf tubes	Multiple Vendors
PCR Grade Water	Multiple Vendors
Qubit™ Fluorometer including High Sensitivity consumables	Thermo Fisher
Thermocycler with heated lid (0.2 mL standard PCR tubes) compatible with semi-skirted ABI style PCR plates and option for ramp rate programming of 0.1 °C/sec (or 2% of max).	Multiple Vendors
The next five items are only necessary when sequencing in-house (not needed in case of sending samples to a core facility)	
NaOH solution (2 N)	Multiple Vendors
Tris/HCl (200 mM), pH 7	Multiple Vendors
Low TE (10 mM Tris-HCl (pH 8.0), 0.1 mM EDTA)	Multiple Vendors
Illumina® NGS Sequencing instrument*	Illumina®
Illumina® Sequencing Reagent kit (300 (2x151) bp run)	Illumina®

\*Note: When running the libraries on a MiSeq System, Exon 11 will not be able to be sequenced bi-directionally due to a sequence motif present in that area of the genome which is incompatible with the MiSeq Chemistry.

## Procedure

### General precaution

Use a Pre-PCR environment for setting up the RC-PCR. Pooling, cleaning and library preparation should be performed in a Post-PCR environment.

### 1. Thermocycling program

Temp	Duration	Ramping rate (from previous step)	Cycles
98 °C	2 minutes	N/A	1 x
98 °C	10 seconds	Max	1 x
58 °C	10 minutes	<b>0.1 °C/sec (or 2% of max)</b>	
72 °C	1 minute	Max	
95 °C	10 seconds	Max	2 x
80 °C	1 second	Max	
62 °C	90 <u>minutes</u>	<b>0.1 °C/sec (or 2% of max)</b>	
72 °C	30 seconds	Max	
95 °C	10 seconds	Max	34 x
80 °C	1 second	Max	
62 °C	2 minutes	<b>0.5 °C/sec (or 10% of max)</b>	
72 °C	30 seconds	Max	

Heated lid at 105 °C.

Double Check: Depending on the instrument, this protocol takes 6-7 hours to complete.

## 2. Reverse Complement PCR

The target specific RC-probes will be transformed into the functional, tailed and indexed PCR primers, followed by multiplex DNA amplification of the target regions.

### 2.1 Thaw on ice:

- RC-PCR Probe Panel A (Black cap)
- RC-PCR Probe Panel B (Red cap)
- Probe Dilution Buffer (Blue cap)
- 2x HiFi Master Mix (White cap)

**Note: The HiFi Master Mix contains isostabilizers and may not freeze completely, even when stored at -15 °C to -25 °C. The Master Mix may contain precipitates when thawed at +2 °C to +8 °C. Always ensure that the Master Mix is fully thawed and thoroughly mixed before use.**

### 2.2. Take two identical IDX PCR plates and cut off the number of strips needed. Mark the plates with 'A' and 'B'.

**Note: Register the indexes used (IDX set/strip-column number and well position for each sample).**

**Note: For each sample, two PCR reactions are needed (Panel A and Panel B). Always use the same well position for the same sample, in order to generate identical indexes for each sample in both panels.**

**Note: Make sure to use the strips or plates in the right orientation.**

### 2.3. Prepare in a fresh 1.5 mL Eppendorf tube the Probe-Polymerase premix Panel A, by combining and mixing:

- 0.2 µL RC-PCR Probe Panel A per reaction (Black cap)
- 1.8 µL Probe Dilution Buffer per reaction (Blue cap)
- 4 µL Molecular Grade Water
- 10 µL Master Mix per reaction (White cap)

### 2.4. Prepare in a fresh 1.5 mL Eppendorf tube the Probe-Polymerase premix Panel B, by combining and mixing:

- 0.2 µL RC-PCR Probe Panel B per reaction (Red cap)
- 1.8 µL Probe Dilution Buffer per reaction (Blue cap)
- 4 µL Molecular Grade Water
- 10 µL Master Mix per reaction (White cap)

**Example: 24 samples + 10% extra volume\***

- **Probe-Polymerase premix:**
  - 5.28 µL RC-PCR Probe Panel
  - 47.52 µL Probe Dilution Buffer
  - 105.6 µL Molecular Grade Water
  - 264 µL Master Mix

**\* It is recommended to allow for a 10% excess when preparing the Probe-Polymerase premix to correct for any pipetting loss. The kit contains extra reagent for this.**

- 2.5. Add to each tube of plate A: 16  $\mu$ L of Probe-Polymerase premix Panel A (from 2.3).
- 2.6. Add to each tube of plate B: 16  $\mu$ L of Probe-Polymerase premix Panel B (from 2.4).
- 2.7. Add to each well 4  $\mu$ L of DNA (5 ng/ $\mu$ L).
- 2.8. Close the tube strips carefully with caps and mix by flicking. Check for a homogeneous pink coloured reaction mix and then spin shortly.
- 2.9. Start the RC-PCR program in the thermal cycler(s) and place the samples in the cycler when the block is between 60 °C and 98 °C. Then close the lid.



**Safe stopping point after RC-PCR**



### 3. Pool, Purify and Sequence

After the PCR, the samples have been amplified and tagged with sample specific indexes and sequencing tails. From this point, PCR products can be pooled together in a single tube, purified by a bead purification to remove primers, primer-dimers and salt.

**Note:** Before pooling, optionally check 3  $\mu\text{L}$  of the unpurified PCR products on agarose (2%).

- 3.1. Bring the AMPureXP™ beads solution to Room Temperature.
- 3.2. Perform steps 3.3 to 3.7 for both Panel A and Panel B individually.
- 3.3. Combine 5  $\mu\text{L}$  RC-PCR products from all the reaction wells from Panel A/B (except negative controls) in a 1.5 mL Eppendorf tube.
- 3.4. Mix well and transfer 40  $\mu\text{L}$  of the pool to a new 1.5 mL Eppendorf tube.
- 3.5. Add 60  $\mu\text{L}$  Low TE buffer or molecular grade water to the tube (total volume is now 100  $\mu\text{L}$ ).
- 3.6. Beads purification:
  - a. Vortex the beads thoroughly to resuspend.
  - b. Add 90  $\mu\text{L}$  beads solution to the 100  $\mu\text{L}$  pool (from step 3.5) and mix well immediately by pipetting up and down 5 times.
  - c. Incubate for 5 minutes, off magnet.
  - d. Place the tube on magnet for 3 minutes or for the solution to be fully cleared.
  - e. Remove and discard all liquid carefully without disturbing the beads.
  - f. Add 200  $\mu\text{L}$  (freshly prepared) 75% ethanol, without disturbing the beads.
  - g. Wait for 1 minute.
  - h. Repeat steps e., f. and g. for a second ethanol wash step.
  - i. Carefully remove all liquid without leaving traces of ethanol. (Optionally a quick spin can be performed, then place the tube back on magnet and remove excess ethanol)
  - j. Dry with open cap for 2-3 minutes at Room Temperature.  
**Do not over-dry.**
- 3.7. Elution:
  - a. On Magnet: Add 120  $\mu\text{L}$  Low TE buffer to the tube.
  - b. Off Magnet: Re-suspend the beads by flicking or short vortexing.
  - c. Incubate for 2 minutes, off magnet.
  - d. Put the tube on magnet and wait for 1-3 minutes or for the solution to be fully cleared.
  - e. Carefully bring 110  $\mu\text{L}$  of the clear solution to a new 1.5 mL Eppendorf tube, making sure not to transfer any of the beads.


2x Ethanol wash


The libraries are now ready for quantification and qualification.

- 3.8. Determine the final concentration of the libraries by a double Qubit (HS) measurement according to manufacturer's manual.
- 3.9. Verify the libraries on TapeStation or Bioanalyzer, according to the manufacturer's protocol. If needed, dilute the panels. Example: For TapeStation High Sensitivity kit, dilute to ~2 ng/ $\mu$ L.

*Examples of clean libraries on TapeStation:*

Panel A





- 3.10. Perform sequencing on an Illumina® NGS platform, according to the manufacturer's manual.
- a. Read at a minimum of 151-8-8-151.
  - b. The Bed. File is available at:  
[https://nimagen.com/gfx/Human\\_Genetics/TP53\\_customer.txt](https://nimagen.com/gfx/Human_Genetics/TP53_customer.txt)
  - c. A minimum read depth per sample of 70 000 will result in a lower detection limit of 5% variant allele frequency.

**Note:** For technical assistance contact our technical support at [techsupport@nimagen.com](mailto:techsupport@nimagen.com).

Apply local regulations regarding waste disposal.

*Legal Notices:*

RC-PCR is patent protected (PCT/GB2016/050558,WO2016146968A1) and exclusively licensed to NimaGen B.V. Nijmegen

Qubit is a trademark of Thermo Fisher Scientific Inc.

AMPureXP is a trademark of Beckman Coulter

**Published by:**

NimaGen B.V.  
Hogelandseweg 88  
6545AB Nijmegen  
The Netherlands  
[www.nimagen.com](http://www.nimagen.com)  
© 2022 NimaGen  
All rights reserved.