

Ready to level up your 16S analysis?

Our EasySeq™ Full-length 16S Library Prep Kit has shown excellent composition representation, high concordance with routine diagnostics – and often detects more species than other 16S-based routine molecular diagnostics. That's why leading microbiology labs are switching to RC-PCR.

“That's why leading microbiology labs are switching to RC-PCR”

Try it yourself!



- ✓ **Sample-to-result within 8 hours**
Ideal for rapid diagnostics and clinical workflows.
- ✓ **Closed-tube, single reaction RC-PCR workflow**
Minimizes pipetting errors and contamination; increases data reliability.
- ✓ **High sensitivity**
Perfect for rare or challenging samples.
- ✓ **High multiplexing: 96 barcodes available**
Increasing throughput - Saving cost - avoids library carryover issues.
- ✓ **Automation-ready & compatible with ONT LSK kits**
Ready for scalable workflows.



Request your sample kit today!

Please contact us:
info@nimagen.com
+31 (0) 24 820 02 41

NimaGen.

NimaGen B.V.

Hogelandseweg 88
6545 AB Nijmegen
The Netherlands

T +31 (0)24 820 0241
E info@nimagen.com
www.nimagen.com

Version 1.1 - April 2025

Evaluation of the EasySeq Full-Length 16S sequencing kit for bacterial identification in clinical samples: a multicentre study

T.T. Severs^{1,2}, J.H.B. van de Bovenkamp³, L.B. van Alphen⁴, L. Gard², J.W.A. Rossen^{1,2,5}, R.Y. Yahiaoui⁶, S.P.F. Matamoros⁷, J. Dingemans⁸, S.D. Pas^{9*}, S.A. Boers^{10*}

¹Laboratory of Medical Microbiology and Infectious Diseases & Isala Academy, Isala Hospital, Zwolle, The Netherlands; ²Department of Medical Microbiology and Infection Prevention, University Medical Center Groningen, University of Groningen, Groningen, The Netherlands; ³Laboratory for Medical Microbiology and Public Health, Labmicta, Hengelo, The Netherlands; ⁴Department of Medical Microbiology, Care and Public Health Research Institute (CAPRI), Maastricht University Medical Center+ (MUMC+), Maastricht, The Netherlands; ⁵Department of Pathology, University of Utah School of Medicine, Salt Lake City, United States of America; ⁶Department of Medical Microbiology, Laboratory, Molecular Diagnostics Unit, Maasstad Hospital, Rotterdam, The Netherlands; ⁷Amsterdam UMC location University of Amsterdam, Medical Microbiology and Infection Control, Amsterdam, The Netherlands; ⁸Laboratory for Molecular Diagnostics, Department of Clinical Biology, Jessa Ziekenhuis, Hasselt, Belgium; ⁹Department of Medical Microbiology, Radboud University Medical Center, Nijmegen, The Netherlands; ¹⁰Center of Infectious Diseases (LUCID), Leiden University Medical Center (LUMC), Leiden, The Netherlands; *Equal contribution

Introduction

Bacterial culture remains an essential method for **pathogen identification** but is often inadequate for fastidious bacteria or complex infections. Its **sensitivity** is further diminished by prior **antibiotic exposure**.

The **molecular** diagnostic revolution provided numerous options for **pathogen detection** or **identification** like syndromic testing by qPCR, Sanger Sequencing or Next-Generation sequencing. Each of these methods offer a different detection spectrum.

In this **consortium**, we compare different molecular methods of 9 different institutes with the **NimaGen's EasySeq Full-Length 16S** sequencing kit in combination with the **Emweb's Genome Detective 16S** pipeline, aiming to assess its suitability for **clinical diagnostics**.

Methods

Nine major medical centres in the Netherlands and Belgium processed **20 clinical samples** each from sterile sites (synovial fluid, cerebrospinal fluid, and tissue), ranging from negative (n=5) to polymicrobial infections.

Each centre applied **16S-based routine molecular diagnostics (MDx)** and the **EasySeq™ Full-Length 16S Library Prep Kit (NimaGen)** for Oxford Nanopore Technologies (ONT) sequencing. All centres worked according to manufacturers protocol and used the same lotnumbers of reagents and MinION flowcells.

Sequence data were analyzed using the Default Bacterial Analysis (16S) pipeline available at **GenomeDetective.com** (Emweb).

To assess **interlaboratory performance** and **reproducibility**, all centres processed Microbial Community DNA Standard II (Zymo, cat# D6310) in triplicate.

Table 1: Contributing centres and their routine 16S-based MDx

Centres	Routine 16S-based MDx
Isala	16S Sanger sequencing
UMCG	16S Sanger sequencing
Labmicta	qPCR
AUMC	IS-pro
LUMC	16S micPCR/Illumina sequencing
RUMC	16S RC-PCR NimaGen, Illumina sequencing
Maasstad ziekenhuis	16S ONT sequencing
Jessa ziekenhuis	16S modified ONT sequencing
MUMC	16S Sanger sequencing

Results




Figure 1: Microbial community DNA standard II performance for all centres tested in triplicate. Centres are indicated by A to I. Centre D only tested in duplo for which was normalized (SEM). 10K reads were subsetting prior analysis to normalize comparison. T indicates the theoretical input of the standard. The blue line in the figure represents 0.01% abundance. With our methods, detecting bacteria with an abundance below this threshold is not feasible.




Figure 2: VENN diagram of bacterial taxa found. 16S-based routine MDx methods are shown in Table 1. Negative Processing Controls revealed bacterial reads (kitome), requiring correction before data visualization. Without this adjustment, taxa uniquely detected by the NimaGen Full-Length 16S kit in combination with the Genome Detective 16S pipeline, would exceed several hundred. Results are shown for 161/180 samples (89.44%) after excluding 19 samples due to failures in library preparation, sequencing and/or data analysis.

Discussion/Conclusion

Based on qualitative results, the EasySeq Full-Length 16S kit shows **high concordance** with routine diagnostic methods. It frequently detects **more bacterial** species than other methods and demonstrates the capability to identify a **high diversity** of bacteria within a single primary sample. The performance of the **Microbial Community DNA Standard** varies slightly between centres but remains consistent within triplicates. Detected bacteria follow the expected order of decreasing absolute load (Figure 1). The **NimaGen Full-Length 16S kit for ONT** in combination with the **Emweb Genome Detective 16S pipeline** may be of great diagnostic value, due to the fast work-up and detection of (multiple) fastidious, clinically relevant bacteria; however, users should be **aware** of the potential influence of **kitome** and **intersample signals**.

More info: t.severs@isala.nl
or [linkedin.com/in/severstt/](https://www.linkedin.com/in/severstt/)



isala



Radboudumc
UNIVERSITY MEDICAL CENTER

