

Product Change Notification

Notification date: 2023-11-14

Implementation date of production: 2024-02-14

This letter is to inform you, of an upcoming product change

Products involved: IDX096-U02(D) 96-well Index plates, compatible with Reverse-Complement PCR Next Generation Sequencing Library Prep kits, for Illumina.

Reason for Change: Performance optimization of individual indices.

Specification of Change: The following indices were replaced:
- IDX096-U02(D) position H01

The new indices for these positions are listed in the table below:

Plate location	Index combination name	Index Name	Oligo Sequence 5' --> 3'	Samplesheet sequence	Index Name	Oligo Sequence 5' --> 3'	Samplesheet iSeq/MiniSeq/NextSeq/HiSeq 3000/4000	Samplesheet sequence MiSeq/HiSeq2000-2500/Novaseq
H01	RPU0104	RP-8700	TCCGGATTCA	TGAATCCGGA	RP-6610	CATACCTATA	TATAGGTATG	CATACCTATA

Change Recognition: The new version IDX-plates can be recognized by the following label on the top of the box. This label of version change will be visible till 3 months after implementation of the new version.

New version IDX-plate!

U02v4

Download new IDX file here:
www.nimagen.com/downloads

The printed text on the plastic plate will also display the version of the IDX-plate.

IDX-files Update Version: The latest version of the EasySeq™ Index Primer Plates, Unique Dual Indexes (UDI) can be downloaded from www.nimagen.com/downloads. In addition to the new version (IDX-U02v4), all previous versions will remain available for download.

